

Discover the Terason Advantage

Exceptional image quality

Quick and accurate visualization of musculoskeletal structures and needle guidance



Ease of use

Reduce learning curve, increase efficiency and patient throughput




Portable and durable


Perform exams in any location

Innovative

Ultrasound technology combined with the versatility of personal computing



Achilles Tendon Insertion



Common Extensor Tendon



Vision • Value • Versatility


Visualize, Diagnose and Treat with Confidence

The Complete Ultrasound Solution

Vision

Terason's vision encompasses a commitment to deliver crystal clear image quality, giving you the ability to easily scan, quickly diagnose, and accurately inject.

TeraVision™ and OneButton™ Optimization features are a culmination of that commitment. The MSK t3000 series provides high quality, patient specific images with sharp contrast resolution and enhanced clarity of muscles, tendons, ligaments, joints and soft tissue.



Carpal Tunnel – Median Nerve

Value

Terason's value means your affordable investment in a Terason ultrasound system is protected for years to come. With multiyear warranties, easy upgrades, trade-ins, technical support team and on-site training, this system will quickly prove itself to be a hassle free, turn key solution to your practice.

Versatility

Terason's versatility in portable ultrasound allows imaging to be easily integrated into your everyday practice. The MSK series of imaging solutions provides utility and functionality that extend beyond point of care diagnosis. It's ready when and where you need it.

- Integrated CD/DVD
- Network and Wireless
- DICOM 3.0 Compliant
- Familiar Windows® Interface



12L5: 5 to 12Mhz.
Linear Array



7L3: 3 to 7 Mhz.
Linear Array



5C2: 2 to 5 Mhz.
Curved Array



12HL7: 7 to 12 Mhz.
Linear (Hockey Stick) Array



Available in a cart or
standalone configuration



Terason, a division of Teratech Corporation
77 Terrace Hall Avenue • Burlington MA 01803
1-866-TERASON • 781-270-4143 • www.terason.com
©2009 Teratech Corporation. All rights reserved. 16-3175 Rev 1

Vision • Value • Versatility